

Clenny®A Family Smart


Congratulations on choosing a Chiesi Italia S.p.A. device. We recommend that you read the instructions for use carefully and keep them for future reference.

Clenny A Family Smart is a device operating on Mesh technology intended for the nebulization of aerosol drugs to support the treatment of the respiratory tract. The operating principle it is based on combines in a single product the best features of compressor technology (rapid and effective nebulization of any type of drug) and ultrasonic technology (absolute silence and portability even during use), making it ideal for any type of use on adults and children, both for in-home care and for emergencies, during travel or outdoors (always respecting the operating conditions). Device not intended for life support or to provide patient monitoring functions. Clenny A Family Smart delivers high performance in a compact device. The special design of the medication chamber promotes efficient drug delivery, reducing the residual volume of medication that remains in the chamber and is not neutralised to less than 0.1 mL. The device can also operate with the medication chamber tilted from the vertical position, provided the medicinal product remains in contact with the mesh.

Clenny A Family Smart is a device intended to be used at home indicated for supporting the aerosol administration of medications for the treatment of respiratory tract diseases. The device is suitable for use in the treatment of paediatric and adult patients (patients without cognitive impairments). Medical devices intended for home use can be used without the presence of healthcare professionals. Use of the device by children and people who are not self-sufficient must always be carried out under the close supervision of an adult in full possession of their mental faculties. No special training is required; however, the user manual contains specific points for the correct use of the device.

DESCRIPTION OF COMPONENTS / CONTENTS

- Medicine Chamber – Vibrating Mesh (silver disc visible inside the chamber)
- Mask with elastic band (1 for adults and 1 for children)
- USB-C to Lightning adapter
- Main unit
- ON/OFF button
- ON/OFF LED light
- Power cable for medicine chamber
- External power connection cable (USB Type C connector)
- Instruction Manuals



ATTENTION: The presence of small components poses a risk of suffocation if swallowed by children, infants or pets. The main unit may get damaged or cause injury. It is therefore recommended to store the product and its components in places where children and pets do not have access. If a child swallows small pieces, consult a doctor immediately.

ATTENTION: BEFORE USING THE DEVICE, CONSULT THE USER MANUAL. WHEN ADMINISTERING THE MEDICATION, ALWAYS FOLLOW YOUR DOCTOR'S INDICATIONS. NEVER DISASSEMBLE THE DEVICE. FOR ANY INTERVENTION, CONTACT THE MANUFACTURER'S AND/OR DISTRIBUTOR'S TECHNICAL SERVICE.

GENERAL WARNINGS

- When opening the packaging, check the integrity of the appliance and accessories, paying particular attention to the presence of damage to the plastic parts, which may make live internal parts accessible, and to breaks and/or peeling of the power cable. In such cases, do not connect the device to any power source. Always carry out these checks before each use.
- Before connecting the appliance, always check that the electrical data indicated on the power source label matches the specific rating plate data (output 5V 1A).
- Follow the safety regulations for electrical equipment and namely:
 - Exclusively use original Accessories and Spare parts supplied by the Manufacturer. Only use original parts of the type indicated in this manual.
 - DO NOT THROW INTO FIRE OR PLACE NEAR FLAMES OR HEAT SOURCES.
 - Never immerse the appliance and its cord in water, as this could cause electrical leaks and shocks during operation and the appliance could cause problems/damage and/or become unusable.
 - Put the appliance on flat and stable surfaces in order to avoid accidental falls or overheating the device or external power source.
 - Do not use the device in environments where flammable anaesthetic mixtures with air, oxygen or nitrous oxide are present.
 - Do not use the appliance in anaesthesia and pulmonary ventilation systems.
 - Do not touch the appliance with wet hands and always prevent the device from coming into contact with liquids.
 - Do not leave the device on and/or connected to the power socket when it is not used.
 - Do not pull the power cord of the main unit to disconnect it from the power source or the medicine chamber, but after turning it off, grasp the device to disconnect it.
 - Store and use the appliance and all its parts in places protected against the weather and away from any sources of heat. After each use, always carry out cleaning (and any disinfection) procedures then store the device away from dust and sunlight.
 - In general, the shared use of accessories, and in particular the medicine chamber, by multiple people is not recommended. We therefore recommend personal use of this device.
- For repairs, only contact the Manufacturer's technical support service or a technical support centre authorised by the manufacturer and request the use of original spare parts. Failure to comply with the above can compromise the safety of the device.
- This device must be used exclusively for its intended use and as described in this instruction manual. Therefore, it must be used as an aerosol therapy system. Any use differing from the intended use is considered to be improper and therefore dangerous and neither the Manufacturer nor the Distributor can be held responsible for damages caused by improper, incorrect and/or unreasonable use or by the connection to electrical systems not compliant with current safety standards.
- The medical device requires special precautions regarding electromagnetic compatibility and must be installed and used in accordance with the information provided with the accompanying documents.
- Children and dependent people must always use the medical device under the strict supervision of an adult in full possession of mental faculties.
- Do not leave the device unattended in places accessible to children and/or persons not in full possession of mental faculties as they may strangle themselves with the power cable.
- The medical device may come into contact with the patient via the nebuliser/masks/

mouthpiece, components compliant with the requirements of regulation ISO 10993-1; therefore, allergic reactions and skin irritations are unlikely to occur. The product and its parts are biocompatible in accordance with the requirements of regulation EN 60601-1.

- The device is very simple to operate and therefore no further explanations are required other than those reported in the following instruction manual.
- ATTENTION:** Do not modify this device without the manufacturer's authorisation. None of the electrical and/or mechanical parts of the device have been designed to be repaired by the user. Failure to comply with the above can compromise the safety of the device and, therefore, of its use.
- Use of the device in environmental conditions other than those specified in this instruction manual may jeopardise the safety and performance of the device. In case of transport/storage at temperatures below 5°C, leave the device in an environment with suitable working conditions before using it.
- The materials in contact with the drugs are thermoplastic polymers with high stability and chemical resistance, widely used in the medical field. However, given the variety and continuous evolution of the medications used, it is not possible to exclude interactions and, therefore, it is suggested to:
 - Always consume the medication as quickly as possible after its opening.
 - Always avoid prolonged contact of the medication with the appropriate tray container (Mesh Chamber) and immediately carry out the cleaning procedures after each application.
 - Should abnormal situations (e.g. vitrification or cracks) in the medication chamber arise, do not introduce any solution and do not perform inhalation. Contact the technical service specifying the method of use and the type of medication used.

- Remember to:
 - Do not use the product to nebulise liquids/substances that are not drugs prescribed by your doctor.
 - Carry out the aerosol treatment by respecting the proper way to perform it, the doses and combinations indicated by the medical specialist, using only the accessory recommended by them, based on the disease to be treated.
- Any serious incidents occurring in connection with the product must be reported to the manufacturer and the competent authority in your country. If you do not have this authority's contact information, please contact the manufacturer (or distributor) whose contact information is given in this instruction manual.

- Do not drop or subject the main unit and/or the medication chamber to strong vibrations. They could be irreversibly damaged or create an electric shock hazard.
- The Manufacturer/Distributor cannot be held liable for accidental or indirect damage, if changes have been made to the device, unauthorised repairs and/or technical interventions have been performed or any of its parts have been damaged by accident, improper use and/or abuse. Any intervention, however minor, on the device immediately invalidates the warranty, and in any case does not guarantee compliance with the technical and safety requirements of Regulation (EU) 2017/745 and the relevant reference standards.
- Do not clean the Mesh with cotton balls, brushes or similar! Never touch the Mesh! Please only follow the cleaning instructions given in the manual.

- Do not clean the Mesh with cotton balls, brushes or similar! Never touch the Mesh! Please only follow the cleaning instructions given in the manual.
- Do not clean the Mesh with cotton balls, brushes or similar! Never touch the Mesh! Please only follow the cleaning instructions given in the manual.
- Do not clean the Mesh with cotton balls, brushes or similar! Never touch the Mesh! Please only follow the cleaning instructions given in the manual.

CONTRAINDICATIONS AND ADVERSE REACTIONS

- Patients allergic to the nebulized drug.
- The device should be used with caution or under the supervision of a physician for patients who may experience breathing difficulties, apnoea, or persistent asthma during the nebulization process.
- The device is not intended for life support or to provide patient monitoring functions.

PREPARATION AND USE

When using it for the first time and always after a long period of non-use, we recommend disinfecting the nebulization chamber and mask and mouthpiece accessories beforehand to avoid bacterial propagation and/or contamination. Then follow the instructions in the Cleaning and Disinfection paragraph of this instruction manual and proceed with the operations described below.

- Do not use the device to inhale water and/or distilled water unless specifically indicated and/or prescribed for medical use, as this could cause symptoms to worsen. Do not drop or apply strong pressure or impact to the medicine chamber to avoid damage that could render it unusable or cause electric shock.
- Do not INSERT AND NEBULISE any substances into/with this appliance other than the drugs and medical devices specifically indicated for aerosol nebulization or which have a specific limitation to the technology that it uses (Mesh).
- Do not USE WITH Medical Devices and/or solutions containing pure or high percentage dissolved hyaluronic acid, essential oils, scented fragrances, etc., or loaded with sediment (e.g. Pure Thermal Waters) not specifically indicated for use with these devices and which may otherwise block or irreparably damage the medicine chamber. When using pure hypertonic solutions (e.g. sea water), perform a thorough cleaning operation after each use to avoid blocking the device.

1 - MEDICATION INSERTION

- Open the medication chamber cover by releasing it with your fingers.
- Insert the prescribed medication(s) in compliance with the instructions of your doctor and the medication manufacturer.
- ATTENTION:** The medicine chamber is equipped with a graduated scale that indicates the maximum volume of drug that can be inserted (word MAX), or 8mL (8cc). It is recommended to never exceed the maximum filling value indicated.
- Close the medicine chamber lid. Do not spill medication/liquid onto the contacts at the base of the medication chamber or the main unit cable. If this happens, dry them and remove them immediately with gauze. Using the device wet can cause malfunction/damage or an electric shock.

2 - PRODUCT PREPARATION

Once the medication is inserted in the chamber, proceed with the following operations:

- Connect the medicine chamber to the power cable of the main unit, paying attention to the "polarized" insertion clearly visible on the connectors of both. Securely connect the two connectors together.

- Then connect the accessory intended for the specific treatment by your doctor to the medicine chamber (Optional Mask or Mouthpiece). Make sure you use the correct size mask based on the patient (L for adults and S for children), in order to ensure maximum comfort and effectiveness of the treatment.



- Then connect the main unit to the chosen power source via the cable with USB Type C connector, taking care to respect the power specifications, as better indicated in the paragraph POWER SUPPLY OF THE DEVICE.

3 - POWER SUPPLY OF THE DEVICE

- The device can be connected to any external power source as long as it has a USB Type C connection and is limited to a nominal output of 5V 1A; therefore, it can be powered by a Power Bank, by the internal batteries of cell phones/tablets/PCs or by regularly certified AC/DC power supplies with a USB connector (with Medical Grade certification). The main unit itself is integrated with a specific USB Type C connector with 5V 1A input, so to power the device you simply need to make the USB connection to a suitable power source. If you use your mobile phone/tablet as a power supply, please connect it directly to the device via the USB Type C adapter, while for older iOS devices, please use the supplied USB Type C – Lightning adapter (operation with devices older than iPhone 11/iPad Gen by is not guaranteed). The device's power consumption during use is minimal and does not interfere with the proper functioning of the device or significantly reduce its battery life. For correct operation depending on the model and brand of the phone/tablet used, it may be necessary to manage the device's USB power settings to enable powering the device (OTG functions and similar).
- None of the power sources mentioned are supplied with the device (power bank or power supply unit). The device meets the electrical grade and safety requirements of IEC/EN 60601-1.

- Regardless of the type of source used, it is important to remember that it is always necessary to ensure that it cannot be damaged by the nebulization of the drug generated by this device. Direct spilling of liquid, or close, direct exposure to a stream of mist, on any electrical device can cause harmful condensation or deposits. Therefore, avoid exposing electronic/electrical devices to the risk of accidental spillage or to the direct and close flow of the mist.

- Operation with AC/DC power supply**
It is important for the safety of the device and its user to only use power supplies of the indicated type (with USB Type C power connector, with 5V 1A output) and ensure that they are certified in accordance with the above.

- To avoid potentially dangerous conditions, overheating, falls and/or accidental damage, it is recommended to always unplug the device from the relative power source immediately after finishing the treatment.

- Always carefully disconnect the device from the pre-selected power supply source and from the nebulization chamber to avoid breakages that compromise its functioning. Any damages resulting from liquid penetration or mechanical failure of the connector are not covered under the Warranty.

- DO NOT FORGET TO TURN OFF THE PRODUCT once you have finished using it. The device DOES NOT AUTOMATICALLY TURN OFF once the medicine runs out, but only after 20 minutes of continuous operation.

4 - USE OF THE DEVICE

Once properly prepared, the device is then ready to be connected to the chosen power source and used as follows:



- During treatment**
 - Breathe calmly and deeply to best inhale the nebulized product
 - For solutions with high surface activity (bubbles) or high viscosity, the nebulization speed may normally decrease.
 - It is normal to hear a slight sound produced by the Mesh during nebulization and until it turns off.
 - If the solution has a high surface activity, it is normal for a small amount of solution to evade from the surface of the Mesh.
 - When the liquid in the chamber runs out, remember to ALWAYS turn off the device to avoid potential damage to the Mesh.
 - When the medication is running out, users are advised to keep the nebulization chamber as vertical as possible, to facilitate the nebulization of the residual liquid by the Mesh.
 - When using the mask accessory, do not cover the mask openings with hands or other objects.
 - When operating the device, avoid blocking/covering the vent hole in the medicine chamber lid, so as not to affect correct nebulization.
 - While inhaling, if the user feels unwell, turn off the device and consult a doctor immediately.

IMPORTANT: The device can be used continuously, but it is important to remember that it is equipped with a protection function against involuntary switching on or failure to switch off, a function that automatically switches it off after 20 minutes of continuous operation (without switching off).

- ATTENTION:** During use, the device must be held vertically or can be tilted in all directions with a wide angle of inclination thanks to the particular conical structure of the medicine chamber, but always enough to allow the drug to reach the nebulization Mesh. If the nebuliser cup is tilted too

much the medication cannot reach the Mesh, the flow of nebulized medication may be interrupted.

- To prevent clogging and/or damaging the Mesh, the product must be operated until the medication is fully used up (until steam output has completely stopped). Completely empty the chamber if treatment is stopped before its entire contents have been used. In both cases, rinse only the medication chamber well under lukewarm running water (low pressure jet to avoid damage to the Mesh component) and insert a few drops of demineralised water. Engage the chamber to the device and switch the device on for the time required (a few seconds) to drain the last residues. When no more steam is output, switch the product off and clean as described below.**

CLEANING AFTER EACH USE

This product must be cleaned after each use, not only for the proper maintenance of the hygienic condition, but to maintain performance and proper operation.

- The failed or incorrect cleaning of the medication chamber immediately after its every use, may partially or completely jeopardise its proper operation and consequently require its replacement (not covered by warranty).

AT THE END OF EACH USE, PROCEED AS FOLLOWS:

- After rinsing and emptying the medicine chamber of any residues of the nebulized substance, start the Mesh Self-Cleaning Function as indicated below:

Phase 1: Place a small amount of clean (preferably distilled or demineralised) water in the medication chamber.

Phase 2: With the device switched off, press and hold the ON/OFF button for at least 5sec. The green light on the front will turn on steady to indicate that the self-cleaning mode has been activated.

Phase 3: After 30 seconds, the self-cleaning function ends and the device automatically switches off.

Phase 4: Discard any water remaining in the medication chamber, switch the device back on in normal mode for a few seconds, i.e. until nothing more is nebulised, indicating that there is no more liquid inside the Mesh. Dry the unit thoroughly with a soft cloth without touching the Mesh.

RUN SELF-CLEANING AFTER EACH USE

- Then proceed with the operations described below:
 - After turning off the device, disconnect the main unit from the medicine chamber and the power source, then store it in a dry and safe place. Then proceed with cleaning the accessories.
- Perform manual pre-cleaning until the accessories are visually clean. Soak them in a cleaning solution. Clean the surfaces with a soft-bristled brush.



- For manual cleaning of accessories (medicine chamber, masks and optional mouthpiece) the following solutions/procedures can be used:
 - 0.5%-1% Neodisher MedClean Forte, soak for 10-30 minutes at 30°C.
 - If using other detergents, select soak times, concentrations, and solution temperatures as directed by the detergent manufacturer.

- Wear disposable protective gloves.
 - Place the accessories in a suitable cleaning basin.
 - Prepare a cleaning solution following the detergent manufacturer's instructions.
 - Completely immerse the accessories in the cleaning solution. Make sure there are no air bubbles on the surfaces of the accessories during immersion.
 - Brush the accessories under the cleaning solution with a soft-bristled nylon brush for 30 seconds until all visible residue is removed. Pay particular attention to the empty spaces in the product. Then, wash the inaccessible parts with the cleaning solution using a syringe. Repeat the operation 5 times. The cleaning process takes about 10-30 minutes.
 - Remove the accessories from the cleaning tray and rinse them under a gentle stream of running water for 30 seconds at room temperature.
 - Place the cleaned and rinsed accessories in a suitable clean container with fresh deionized water. Rinse thoroughly for 5 minutes in water.

Rinse the Mesh carefully when using it to nebulize particularly dense drugs and with high viscosity such as sublingual agents or expectorants, to remove deposits that often remain on the Mesh.

The mesh is the silver disc visible inside the medication chamber. NEVER touch the Mesh directly. Be sure to clean and store the device after inhalation. Otherwise, the Mesh may become damaged and no longer usable.

- DO NOT CLEAN THE MAIN UNIT USING WATER. AVOID CONTACT WITH WET SURFACES AND EXPOSURE TO POSSIBLE SPLASHES. DO NOT CLEAN THE DEVICE, ITS PARTS OR ACCESSORIES IN A DISHWASHER. DO NOT USE DETERGENTS, CHEMICALS OR SCALE REMOVERS TO CLEAN THE DEVICE AND ITS ACCESSORIES.**

- Place the medication chamber and accessories on a clean and dry cloth and leave to dry completely in the open (about 2 hours) before storing or re-using them.

- The mesh is the silver disc visible inside the medication chamber. Do not use cotton balls, brushes or other objects to clean or dry the medication chamber Mesh. Never touch the Mesh. If necessary, use a dry cotton cloth to dry the accessories and the medication chamber.**

- To clean the main unit and its cables and connectors, use a soft, dry cloth. If necessary, clean and/or disinfect the main unit (including cables) with a soft cloth moistened with a 75% denatured alcohol solution. Always keep the main unit and its cables and connectors clean and dry and remove any drug residue.

- Once air drying is complete, store everything in the package or travel bag together with the accessories. Keep it all in a clean environment and protected against the weather.

DISINFECTION

Before proceeding with the disinfection operations, always follow the cleaning procedure described in the previous paragraph. The MEDICINE CHAMBER and MOUTHPIECE (optional) can be disinfected by immersion in boiling water or water at a temperature equal to or greater than 95°C for a period of 2 minutes, proceeding as follows:

Phase 1: Immerse the accessories in boiling water and let them boil (over low heat) for 2 minutes.

Phase 2: Once the disinfection time has elapsed, carefully remove the accessories from the boiling water.

The disinfection process at 95°C for 2 minutes has been validated to ensure effective disinfection.

MASKS can be disinfected by immersion in a 75% denatured alcohol solution for 2 minutes or by carefully cleaning the surface with gauze soaked in the same solution.

Phase 3: Pat the surfaces of the accessories dry with a sterile, lint-free cloth until no there is no water (or alcohol solution).

Phase 4: Store treated products in a dry, clean, dust-free environment at moderate temperatures, following the information in the instructions for use and on the label.

- Do not use microwaves to sterilize the device and its accessories.**

MAINTENANCE

The appliance does not contain any parts requiring maintenance and/or lubrication. However, it is necessary to perform a few simple checks to verify the functionality and safety of the appliance before each use. Remove the appliance from the box and always check that there is no visible damage; Pay particular attention to cracks in the plastic that may expose some electrical components. Verify that the medication chamber has not been broken from previous use. If any part is damaged, do not use the device. Using a damaged device may result in improper handling, damage, injury and/or serious hazards.

SPARE PARTS PURCHASABLE SEPARATELY

Medicine chamber	REF SP.0325
Mask KIT (Adult + Child with elastic)	REF SP.0326
Mouthpiece with L-connector	REF SP.0327
USB Type C to Lightning Adapter	REF SP.0328 (Customer Service only)

Only use original accessories/spare parts provided/indicated by the Manufacturer

MEDICATION CHAMBER: For each individual patient, it is recommended to replace the medicine chamber 3 months after the first treatment. In addition, the medication chamber must be replaced in the case of cracks or if the holes in the Vibrating Mesh inside it are partially or totally obstructed by medication, dust, limestone deposits, etc. In the presence of diseases with risk of infection and microbial contamination, we recommend a personal use of the accessories and medication chamber and compliance with the disinfection procedures described in this instruction manual (always consult your medical specialist).

Mask and mouthpiece (the latter is not included in the package) must be replaced if there are visible signs of decay in the material they are made of. Minimum expected service life: 2 years (under standard test and operating conditions).

POSSIBLE PROBLEMS AND SOLUTIONS

Defect	Cause	Remedy
1. Very low or no nebulisation (with device on)	The holes in the vibrating Mesh, located inside the medication chamber, are clogged	Clean the medication chamber as indicated in the instruction manual. If the problem persists even after cleaning, replace the medication chamber with a new one
	The device's power source is running low	Replace and/or restore the power source
	Very oily medicine.	Dilute the medication with a physiological solution
2. No nebulization with the device on (green power LED on)	The electrodes of the medicine chamber or main unit are dirty and/or oxidized	Clean the contacts by removing the cause of the poor electric contact
	Main unit contacts are wet with water or medication	Remove liquid from contacts
	The Mesh is extremely Oxidized or broken	Replace the medicine chamber with a new one
	The main unit cable is faulty	Replace the device
3. No nebulisation and operation LED light off	The main unit is not properly connected to the power source	Check for correct USB connection to a proper power source
	The power source is Completely discharged	Replace and/or restore the power source
	Damaged adapter or disabled external smartphone connection	Replace connector or enable the external connection function on the phone

Defect	Cause	Remedy
4. Fluid leak from the Mesh / medicine Chamber	The medicine chamber is damaged	Replace the medicine chamber
Defects 1 – 2 – 3 – 4	None of the solutions are effective	Contact your dealer or the Manufacturer's technical service

ATTENTION: The Manufacturer declares that the information contained in this instruction manual is consistent with the technical and safety specifications of the device to which the manual refers. The stated technical specifications are updated to the date of the publication of this document and only apply to the machine to which they are attached. The Manufacturer reserves the right to make changes or improvements without prior notice in this document.


TECHNICAL SPECIFICATIONS

TYPE (Regulation 2017/745)	Medical Device with Risk Class IIa
MODEL	Clenny A Family Smart – REF RE 300912/04
POWER SUPPLY	USB Type C input 5V port $\pm 1A$ from external power source
POWER CONSUMPTION	2.0 W approximately
RESIDUAL VOLUME (2 spitting)	< 0.1 mL
MMAD** (ISO 27427:2023)	3.62 (± 0.01) μm
GSD** (ISO 27427:2023)	1.51 (± 0.02)
AEROSOL OUTPUT (ISO 27427:2023)	1.98 (± 0.06) mL
AEROSOL OUTPUT RATE (ISO 27427:2023)	0.21 (± 0.01) mL/min
RESPIRABLE FRACTION (q ₅ μm)	approx. 78.5% (2–5 μm = 70.8±1.39% / 2 μm 57.5±0.70%)
AVERAGE NEBULIZED VOLUME	0.30 mL/min (internal tests with 2 mL NaCl 0.9%)
INSULATION CLASS	Class II (when used with external power source)
WEIGHT	approx. 29 g (without accessories)
SIZE	77 (L) x 14 (H) x 22 (D) mm (total length of cables: 1m)
OPERATION	Continuous
MAXIMUM SOUND LEVEL	± 35 dB (A)
MAXIMUM MEDICINE CAPACITY	8 mL
OPERATING CONDITIONS:	
Room temperature	10 ± 40 °C
Relative ambient humidity	30 ± 85% RH (non-condensing)
Atmospheric pressure	700 ± 1060 hPa
STORAGE AND TRANSPORT CONDITIONS	
Room temperature	-20 ± 60 °C
Relative ambient humidity	30 ± 85% RH (non-condensing)
Atmospheric pressure	700 ± 1060 hPa

*MMAD = Mass Median Aerodynamic Diameter (with AC). **GSD = Geometric Standard Deviation

PARTICLE DISTRIBUTION GRAPH (TEST as per ISO 27427:2023)

The cumulative size distribution graph below shows the results of the drug mass fraction collected at each stage with the Marple Personal Cascade Impactor.



WARRANTY

The product is covered by a 5-year warranty against any defects in conformity with normal conditions of use as set out in the instructions. The warranty does not therefore apply to accidental damage, damage resulting from improper use and wear or consequent accidents. The duration of the warranty against defects of conformity is subject to the specific provisions dictated by national legislation. The medicine chamber and accessories are considered consumables and therefore, if already used, are not covered by the device warranty.

In case of malfunction, contact the TOLL-FREE NUMBER - CUSTOMER SERVICE ITALY 800 005 525 Available on weekdays from 9:00am to 12:00pm and from 1:30pm to 4:30pm (Free when calling from Italy).

Email CUSTOMER SERVICE – assistenza@ca-mi.it

SYMBOLS

Symbol	Definition	Symbol	Definition
	BF type applied part		Attention
	Consult the instruction manual		Serial Number
	Product Identification Code		Production lot
	On / Off		Device with insulation class II
	Alternating current		Manufacturer
	Direct Current		Mark of Conformity to Regulation (EU) 2017/745 on medical devices
	Medical Device		Unambiguous device identification
	Date of manufacture		Distributor
	Humidity limit of use		Temperature limit during use

Symbol	Definition	Symbol	Definition
	Atmospheric pressure limit of use		Temperature limit for storage and transportation
	Humidity limit for storage and transportation		Atmospheric pressure limit for storage and transportation
IP22	Protection rating of an electrical device against accidental or intentional penetration of foreign solid bodies, and protection against the penetration of liquids		
	1st DIGIT: Penetration of solids: Protected against solids having a dimension greater than Ø 12mm		
	2nd DIGIT: Penetration of liquids: Protected against water drops deflected up to 15°		

Device Made in China



WARNINGS REGARDING THE CORRECT DISPOSAL OF THE PRODUCT ACCORDING TO THE EUROPEAN DIRECTIVE 2012/19/EU/WEEE:

After use, the product must not be disposed of with household waste. It can be delivered to specific separate collection centres provided by the municipal administration, or returned to the distributor when purchasing a new device of the same type and having the same functions. Disposing of the product separately avoids possible negative effects on the environment and human health and allows for the recovery of the materials of which it is composed in order to obtain significant energy and resource savings. The symbol on the data label indicates the separate collection of electric and electronic equipment.

Attention: Incorrect disposal of electrical and electronic equipment may result in fines.

EMC – ELECTROMAGNETIC COMPATIBILITY GUIDE AND MANUFACTURER'S DECLARATION

This section contains information on the compliance of the device with EN 60601-1-2:2015 + A1:2021. Clenny A Family Smart is a medical device that is particularly suitable for home applications. CISPR group and category classification: group 1, category B.

- This equipment must be installed and put into service in accordance with the information provided in the ACCOMPANYING DOCUMENTS.
- IMMUNITY TEST LEVELS for basic safety and essential performance of ME equipment and ME systems should be chosen based on a high probability of maintaining basic safety and essential performance, and should be consistent with the professional healthcare facility environment, home healthcare environment, and special environments, based on the intended locations of use.
- HOME HEALTHCARE SETTING is the place of residence where a patient lives or other places where patients are present, excluding settings in professional healthcare facilities where medically trained staff are continuously available when patients are present. This also includes schools, open spaces, homes, hotel vehicles and guesthouses.